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Cat. No.: 10020282

For microbiology

The kit is designed to test antimicrobial susceptibility of Gram positive bacteria enterococci and A,B,C,G streptococci on the basis of bacterial growth at breakpoint concentrations. The kit contains 40 detection strips.

Breakpoint concentrations and interpretation of sensitivities were determined on the basis of the European standard EUCAST (www.eucast.org) dated from 5.1.2011. The test is based on the rehydration of antibiotics in the wells with Suspension medium for SENSILATEST and addition of bacterial suspension. The results are read visually after 18-24 hours of incubation. We recommend to use disk or gradient tests for antibiotic susceptibility testing of viridans streptococci and *Streptococcus pneumoniae*, due to its poor growth in dilution systems.

Please be aware, that inoculum preparation and interpretation differs for enterococci and streptococci!

The kit contains:

- 10 plates for examination
- 1 pc of Suspension medium for SENSILATEST in dehydrated form – Attention: It is necessary to prepare it before testing!
- A lid
- 10 pc of PE bags
- Record sheet
- Storage bag

Storage and expiration of the kit:

It is recommended to store the kit at (+15 to +25) °C. The date of expiration is indicated on each package. After the aluminium package is opened, it is possible to insert any unused strips into a storage bag with desiccant at the room temperature for a period of maximum 8 days. Mark the storage bag with a type of stored strips (e.g. G-I, G-II etc.) to avoid mistake in their further usage. Any longer exposition to open air leads to antibiotic activity failure!!! In case that laboratory is not able to keep kits at recommended room temperature (not exceeding +25 °C), it is better to keep kits in refrigerator (+2 to +8) °C. Then it is necessary to leave plate at room temperature for 2 hours before the use to avoid water condensation.

Storage and expiration of Suspension medium for SENSILATEST:

It can be stored in powder form at (+15 to +25) °C till the date indicated on the package. The expiration of the Suspension medium for SENSILATEST is 6 months at (+15 to +25) °C after it was prepared into a liquid form. Do not store the prepared medium at lower temperatures as it can lead to formation of gel in the media. If such a situation occurs, adjust the medium to room temperature and shake it thoroughly before the use.

Material required to perform a test, not included in the kit:

- Sterile physiological solution
- Sterile tubes
- An stepper pippete for dosage of 100 µl
- A pippete for dosage of 40 µl and 80 µl
- Densitometer (e.g. Densi-La-Meter II, Erba Lachema)
- Incubator 37 °C
- Regular microbiological laboratory equipment (loops, marker, burner, etc.)

Material required to prepare Suspension medium for SENSILATEST, not included in the kit:

- A screw top bottle
- A mixing stick or a stirrer
- Autoclave

Caution: *The kit is for professional use only!
Respect the rules for work with infectious material!*

Instructions for Use

Caution: *Prepare the Suspension medium before you start testing!*

Preparation of Suspension medium:

One bag is sufficient to prepare 0.5 L of Suspension medium for SENSILATEST. Pour the bag content in to a screw top bottle. Add 0.5 L of preheated distilled water and stir it. Autoclave the prepared suspension for 30 minutes at 121 °C. Cool it to room temperature and stir it again under sterile conditions. A small precipitate at the bottom does not influence the function of the medium. Check the pH of the medium before the use. pH of the medium must be 7,5 (± 0,1). Adjust pH if necessary.

Instructions to use the kit:

- 1) Prepare a tube with 2 - 5 ml of physiological solution of pH 5,5 – 6,5.
- 2) Prepare a tube with Suspension medium for SENSILATEST optimally preheated to 37 °C
- 3) Remove few colonies from 18 – 24 hour culture on blood agar and prepare a bacterial suspension in physiological solution of density:

for enterococci:

0,5 on McF scale

for streptococci:

1 on McF scale

Overestimation of suspension density is smaller error than its underestimation.

Inoculation:

Remove a plate from aluminium bag. Remove required number of strips for examination. Do not remove the strips from the package earlier than 30 minutes before you start the work.

Insert any unused strips into a ZIP log bag. Transfer a desiccant bag from the original package into a ZIP log bag and close it for later use. Use the strips stored in this way within 8 days! Mark the storage bag with a type of stored strips (e.g. G-I, G-II etc.) to avoid mistake in their further usage.

Any longer exposition to open air leads to antibiotic activity failure!!!

Record numbers of the examined strains on the corresponding strips. Mark the frame with a type of strip (e.g. G-I, G-II etc.) to avoid mistake in reading results after incubation.

Enterococcus:

Inoculate **40 µl** of bacterial suspension into 4.4 ml of Suspension medium for SENSILATEST and homogenise well. Inoculate 100 µl of suspension in suspension medium into each well of the strip with the help of automatic stepper pipette.

Streptococcus:

Inoculate **80 µl** of bacterial suspension into 4.4 ml of Suspension medium for SENSILATEST and homogenise well. Inoculate 100 µl of suspension in suspension medium into each well of the strip with the help of automatic stepper pipette.

Incubation:

Insert the inoculated plate into a PE bag. Fold the open end of the bag under the plate to prevent evaporation during the incubation.

Incubate the plate at 37 °C for 18 – 24 hours.

Evaluation:

Remove the plate from the PE bag. To read the growth in the microwells, choose a way which is the most convenient for you:

- Use the lid with imprinted antibiotics sequence.
- Read against a grey background.
- Read against natural or artificial dispersed light – you can detect even limited growth in a form of grains or light turbidity.

Please read with attention:

You must see a growth in a control well! If the growth is not present, the test MUST NOT be evaluated! However, the growth control is not used to compare the presence of the growth in individual wells with antibiotics. Every change from a total transparency is evaluated as growth, including small grains of growth in a well, light turbidity or partial growth covering only a part of a well. Beware to differentiate grains of growth from media bubbles. Record the results.

H	G	F	E	D	C	B	A
Growth control	Chloramphenicol	Moxifloxacin		Linezolid		Ampicillin	
GC	CMP 8	MOX 0.5	MOX 1	LIZ 2	LIZ 4	AMP 4	AMP 8
Tigecycline		Levofloxacin		Erythromycin		Penicilin G	
TGC 0.25	TGC 0.5	LEV 1	LEV 2	ERY 0.25	ERY 0.5	PEN 0.25	PEN 2
Ceftriaxon		Tetracycline		Clindamycin	Gentamicin	Vancomycin	
CRO 0.5	CRO 1	TET 1	TET 2	CLI 0.5	GEN 256	VAN 2	VAN 4

Interpretation:

Growth control	GC	Enterococcus										Streptococcus A,B,C,G																	
Ampicillin	AMP	-	-	S	+	-	I	+	+	R	See note below																		
Ceftriaxon	CRO	Don't interpret															-		S				+		R				
Clindamycin	CLI																-	-	S	+	-	I	+	+	R				
Erythromycin	ERY	Don't interpret										-	-	S	+	-	I	+	+	R									
Gentamicin	GEN											-		S				+		R	Don't interpret								
Chloramphenicol	CMP	Don't interpret										-		S				+		R									
Levofloxacin	LEV											-	-	S	+	-	I	+	+	R	-	-	S	+	-	I	+	+	R
Linezolid	LIZ	-	-	S	+	-	S	+	+	R	-	-	S	+	-	I	+	+	R										
Moxifloxacin	MOX	Don't interpret										-	-	S	+	-	I	+	+	R									
Penicilin G	PEN											-	-	S	+	-	R	+	+	R	-	-	S	+	-	I	+	+	R
Tetracycline	TET											-	-	S	+	-	I	+	+	R	-	-	S	+	-	I	+	+	R
Tigecycline	TGC	-	-	S	+	-	I	+	+	R	-	-	S	+	-	I	+	+	R										
Vancomycin	VAN	-	-	S	+	-	S	+	+	R	-	-	S	+	-	R	+	+	R										

If there is a growth at a higher concentration of an antibiotic and no growth at a lower concentrations of the same antibiotic, do not evaluate susceptibility to this antibiotic. This combination shows an error in the procedure !

Enterococci:

- Ampicillin:** result can be used also for AMO, AMO/CLA, AMP/SUL, PIP, PIP/TAZ (these antibiotics are not included in the set)
- Gentamicin:** evaluate in infectious endocarditis only
- Linezolid:** if resistance found confirme by MIC gradient (e.g. E-test). Resistant strains are rare
- Tigecycline:** if resistance found confirme by MIC gradient (e.g. E-test). Resistant strains are rare
- Vancomycin:** if resistance found confirme by MIC bradient (e.g. E-test) or other method (PCR)

Streptococci A,B,C,G:

- Clindamycin:** in case of ERY resistant along with CLI sensitive or intermediate, it is necessary to determinate „induction“ of MLSb (macrolides-lincosamides-streptogramins B) resistance with discs ERY and CLI and correct the result for CLI.
- Erythromycin:** results can be used also for azithromycin, clarithromycin, roxithromycin
- Levofloxacin:** resistance is necessary to be confirmed by disk diffusion test or by gradient MIC (e.g. E-test)
- Linezolid:** resistance is necessary to be confirmed by disk diffusion test or by gradient MIC (e.g. E-test)
- Moxifloxacin:** resistance is necessary to be confirmed by disk diffusion test or by gradient MIC (e.g. E-test)
- Penicilin:** Susceptibility to ampicillin ,cephalosporins (ceftriaxon) and carbapenems (not included in this kit) are inferred from the penicilin. Ampicillin and ceftriaxon can be used however as a control: in case that strain is susceptible to penicilin, results in ampicillin and ceftriaxon should be - -, if not error is indicated. Strains resistant to penicilin are not described so far or are very rare – these strains should be sent to reference laboratory for confirmation. Betahemolytic streptococci don't produce betalactamases – combination with inhibitors thus not provide any benefit.
- Tigecycline:** if resistance found confirme by MIC gradient . Resistant strains are rare
- Vancomycin:** if resistance found confirme by MIC gradient or other method (PCR)

Depending on national or laboratory standards, other criteria described in EUCAST Expert Rules in Antimicrobial Susceptibility Testing can be used.

Quality control:

Control strains can be used for internal testing of functionality of the antibiotics in the laboratory. Please refer to the acceptable limits for quality control strains according to EUCAST.

CCM 4223 <i>Staphylococcus aureus</i> (ATCC 29213)							
GC +	CMP +/- 8	MOX - 0,5	MOX - 1	LIZ +/- 2	LIZ - 4	AMP - ¹ 4	AMP - ¹ 8
TGC - 0,25	TGC - 0,5	LEV - 1	LEV - 2	ERY +/- 0,25	ERY +/- 0,5	PEN +/- 0,25	PEN - 2
CRO ² 0,5	CRO ² 1	TET - 1	TET - 2	CLI - 0,5	GEN - 256	VAN - 2	VAN - 4

CCM 3954 <i>Escherichia coli</i> (ATCC 25922)							
GC +	CMP - 8	MOX - 0,5	MOX - 1	LIZ ² 2	LIZ ² 4	AMP +/- 4	AMP +/- 8
TGC - 0,25	TGC - 0,5	LEV - 1	LEV - 2	ERY ² 0,25	ERY ² 0,5	PEN ² 0,25	PEN ² 2
CRO - 0,5	CRO - 1	TET ² 1	TET ² 2	CLI ² 0,5	GEN - 256	VAN ² 2	VAN ² 4

1 no result stated in EUCAST, result stated by internal quality control Erba Lachema

2 no result stated in EUCAST

+ growth

- no growth

ATCC – American Type Culture Collection
CCM – Czech Collection of Microorganisms

Health protection:

Components of the kit do not contain any dangerous substances.

Disposal of the used material:

Insert the used plate into the vessel intended for the infectious material and autoclave or destroy it by incineration.
Put paper packaging waste to recycling.

Date of last revision: 7.9.2011